


Optimization of HTS Assay Development with Integration of SpyTag/SpyCatcher Technology

MARIAH CULPEPPER, VANESSA CHIRINO, CERISE RAFFIER, CHRISTOPHER LEMKE


Center for the Development of Therapeutics, Broad Institute of MIT and Harvard, Cambridge, MA, USA

ABSTRACT

The Cancer Dependency Map¹ identified protein A as a synthetic lethal target in drug-resistant cancer cells. We aim to identify small molecules that disrupt the A-B-C complex.




A Time-resolved Fluorescence Energy Transfer (TR-FRET) assay was developed to do high-throughput screening (HTS) for hit identification. When protein C is displaced by a small molecule, the signal will decrease.



An ideal labeling strategy would need to bring the donor/acceptor pair close enough to signal, be small enough to minimize background fluorescence, and be cost-effective for HTS. Labeling strategies are used because the donor/acceptor fluorophores cannot bind directly to the proteins.

Past commercially available labeling strategies did not meet these parameters, but **SpyTag/SpyCatcher technology** was made in-house to favorably link protein C and the acceptor (ATTO647) fluorophore. SpyCatcher covalently binds to SpyTag which is linked to protein C, increasing the signal-to-background ratio in the assay.





OBJECTIVES

1. Express, purify and characterize the SpyTag/SpyCatcher technology.
2. Optimize the labeling reaction of protein C-SpyTag-SpyCatcher and the ATTO647 acceptor dye for HTS.


EXPRESSION AND VALIDATION OF SPYTAG/SPYCATCHER

1. How do we characterize the covalent attachment of SpyTag/SpyCatcher to protein C?

SpyTag/SpyCatcher System²




Methods³



Covalent binding allows for more consistent dye labeling and stronger TR-FRET signal

Covalent binding between Protein C-SpyTag-SpyCatcher




SDS-PAGE gel stained with Coomassie Blue

OPTIMIZATION OF ACCEPTOR DYE LABELING


2. How can we optimize the conditions of ATTO647 labeling to be more time and cost-effective in a high-throughput screening format?

Thiol-Maleimide Reaction⁴




Cysteine mutation (SpyCatcher S49C) allows for a thiol-maleimide reaction, labeling ATTO647 on protein C-SpyTag-SpyCatcher

Methods⁵



ATTO647 Labeling of Protein C-SpyTag-SpyCatcher



Acceptor dye use reduced four fold
Time of reaction reduced by a factor of 8

Data interpreted from LC-MS/MS

FUTURE DIRECTIONS


TR-FRET and SpyTag/SpyCatcher technology is effectively utilized in hit identification and optimization.

The next steps in the **cancer drug discovery process** will be applying these techniques to an HTS assay in order to screen 500k compounds.

The TR-FRET assay with its current parameters remains high quality in automation, technology used for high-throughput screening.

This assay optimization will limit the number of molecules screened and reduce false positive hits.

Drug Discovery Process⁶



REFERENCES



ACKNOWLEDGEMENTS

I would like to express my gratitude to Vanessa Chirino and Cerise Raffier for their dedicated mentorship and guidance, as well as the entirety of the CDoT team. It has been an honor conducting drug discovery research with you all.

Funding provided by the Broad Institute and NIH Research Grant R25 HG006682 NHGRI